
HearingAids.md

THE INDEPENDENT CONSUMER GUIDE

Hearing Aids, Explained

How to understand your hearing, compare your options and choose with confidence



FREE DIGITAL EDITION | VERSION 1.0

Created exclusively for HearingAids.md | August 2026

Before you begin

A clear, independent starting point - with important limits.

Education, not medical advice

This guide provides general consumer education. It does not diagnose hearing loss, recommend a specific treatment or replace advice from a physician, ear-nose-throat doctor (ENT), audiologist or other qualified hearing professional. If you have urgent or sudden symptoms, seek prompt medical care.

What this guide will help you do

- ☐ Recognize common signs of hearing difficulty and know when professional care matters.
- ☐ Understand the main differences between over-the-counter (OTC) and prescription hearing aids.
- ☐ Compare features, support, total cost, trial terms and warranties without sales pressure.
- ☐ Prepare for the first month of use and solve common day-to-day problems.

Publisher disclosure

HearingAids.md is an educational consumer resource operated by Sonivue Lab. Product specifications, availability, prices, benefits and policies can change. Verify current information with the manufacturer, seller, insurer and your hearing professional. Product mentions do not guarantee suitability or results. If affiliate links appear, they should be identified where used.

Editorial basis: official FDA, NIH/NIDCD and Medicare sources listed on page 26. Individual results vary.

Contents

Use this guide from beginning to end, or jump directly to the decision you are making.

01 **Types of hearing loss and your starting point** 4-8
Practical explanations, checklists and worksheets

02 **What to expect with hearing aids** 9-13
Practical explanations, checklists and worksheets

03 **Choosing a hearing aid and buying wisely** 14-18
Practical explanations, checklists and worksheets

04 **Living well with hearing aids** 19-23
Practical explanations, checklists and worksheets

05 **Your action plan and references** 24-28
Practical explanations, checklists and worksheets

TIP

Print pages 12, 14, 18, 24 and 27 as worksheets. Bring the completed pages to an appointment or keep them open while comparing products online.

Start with the right path

Your symptoms, age and support needs matter more than any brand name.

1

Check for warning signs

If any red flag on page 7 applies, contact a doctor - preferably an ENT - before buying a hearing aid.

2

Understand your hearing

A professional hearing evaluation gives the clearest picture. For straightforward perceived mild-to-moderate loss, adults may also consider OTC options.

3

Choose a care route

Decide whether you want self-directed OTC access, professional fitting and follow-up, or a blended option with remote support.

4

Compare the whole package

Look beyond device price: support, return terms, warranty, repair, accessories and ease of use all affect value.

5

Test in real life

Use the trial period to evaluate conversations, television, phone calls, noisy places and comfort - not just loudness.

Based on FDA consumer guidance for obtaining and evaluating hearing aids [2][3].

Types of hearing loss

Hearing is more than volume. Clarity, pitch, timing and listening effort all matter.

Conductive hearing loss

Sound is blocked or reduced in the outer or middle ear. Causes can include earwax, fluid or eardrum problems. Some causes may be medically treatable.

Sensorineural hearing loss

Damage involves the inner ear or hearing nerve. Age, noise exposure, genetics, illness and some medicines can contribute. Hearing aids often help, but do not restore normal hearing.

Mixed hearing loss

A combination of conductive and sensorineural loss. Professional assessment helps separate the components and identify next steps.

One-sided or uneven loss

A noticeable difference between ears or tinnitus in only one ear is an FDA warning sign. Seek medical evaluation rather than self-treating.

What people often notice first

- Speech sounds muffled, especially consonants.
- Conversation in restaurants or groups takes extra effort.
- Television volume seems comfortable to you but loud to others.
- You ask people to repeat themselves.
- Phone calls are harder than face-to-face conversation.
- You avoid gatherings because following speech is tiring.

Sources: NIDCD hearing-aid and age-related hearing-loss guidance [5][6].

Is it time for a hearing check?

A simple self-check can reveal patterns, but it cannot diagnose the cause or degree of hearing loss.

- [ ] Do you often increase television or radio volume?
- [ ] Is speech difficult to follow at parties, restaurants or family gatherings?
- [ ] Do you feel frustrated or left out during conversations?
- [ ] Do you avoid groups or social visits because listening is difficult?
- [ ] Do you have trouble hearing on the phone?
- [ ] Does listening leave you unusually tired?
- [ ] Do people seem to mumble even when they are nearby?
- [ ] Has someone close to you expressed concern about your hearing?

What your answers mean

Several 'yes' answers are a good reason to arrange a hearing evaluation. A test can identify the type and degree of loss, measure speech understanding and help determine whether medical care, hearing aids or another strategy makes sense.

Who can help?

A primary care clinician can consider general medical factors. An ENT diagnoses and treats medical or surgical ear conditions. An audiologist evaluates hearing and provides non-medical hearing rehabilitation, fitting and follow-up.

Self-check adapted in plain language from NIDCD consumer guidance [6].

SAFETY FIRST

Stop and seek medical care

Do not rely on an OTC hearing aid alone if any warning sign below applies.

FDA warning signs

- A birth difference, unusual ear shape, injury or deformity.
- Blood, pus or other fluid from the ear within the past six months.
- Ear pain or persistent discomfort.
- A lot of earwax, or concern that an object may be in the ear.
- Severe dizziness, spinning, swaying or repeated vertigo.
- Hearing that changed suddenly, worsened quickly or fluctuated within the past six months.
- Much worse hearing in one ear, or a noticeable difference between ears.
- Ringing or buzzing that occurs only in one ear.
- Difficulty hearing loud sounds or speech even in a quiet room.


Also remember

OTC hearing aids are for adults age 18 and older with perceived mild-to-moderate hearing loss. Anyone younger than 18 needs a prescription pathway and medical evaluation. Sudden hearing loss can be time-sensitive; contact a medical professional promptly.

Source: FDA OTC hearing-aid warnings and benefits/limitations guidance [2][4].

What a hearing test can tell you

An audiogram plots the quietest sounds you can detect at different pitches.



Illustrative only - not a diagnostic audiogram

The axes

Horizontal: pitch or frequency, from low to high.

Vertical: hearing level. Softer sounds are near the top; louder sounds are lower.

The symbols

A common convention uses O for the right ear and X for the left. Your report may include additional symbols for bone conduction and masking.

The bigger picture

Thresholds do not show every real-world difficulty. Speech testing, hearing history, comfort, listening goals and medical findings also matter.

Questions to ask after your test

What type and degree of hearing loss do I have? Is one ear different? How well do I understand speech? Are there signs that require medical care? Would OTC or prescription devices fit my needs? How will benefit be verified after fitting?

Sources: NIDCD hearing-aid guidance [5] and HearingAids.md audiogram education [10].

What to expect with hearing aids

What hearing aids can - and cannot - do, and why realistic expectations matter.



A hearing aid should support communication - not simply make everything louder.

HEARING AIDS MAY

- Make some sounds easier to detect.
- Improve access to speech in many settings.
- Support phone, media and face-to-face communication.
- Reduce listening effort for some people.

HEARING AIDS DO NOT

- Restore normal hearing.
- Remove all background noise.
- Guarantee perfect speech understanding.
- Work equally well for every person or every environment.

Expect an adjustment period

Your voice, dishes, paper and environmental sounds may seem unusually noticeable at first. Regular use, realistic goals and appropriate fine-tuning help the brain adapt. Contact support if sound is painful, distorted or consistently unhelpful.

Source: FDA hearing-aid benefits and limitations [4]. Image is illustrative.

Common hearing-aid styles

Style affects visibility, handling, battery size, features, ear-canal feel and available power.

BTE — Behind-the-ear

Case sits behind the ear and connects to an earmold or tip. Flexible and often easier to handle; available across a wide range of hearing needs.

RIC — Receiver-in-canal

A small case sits behind the ear with a thin wire to a receiver in the canal. Popular, discreet and flexible; receiver and dome need routine care.

ITE — In-the-ear

Custom shell fills part or most of the outer ear. Controls and batteries may be easier to handle than very small devices.

ITC — In-the-canal

Custom device fits partly in the canal. More discreet, but small controls and reduced battery space may affect usability.

CIC — Completely-in-canal

Fits deeper in the canal. Very discreet, but may be harder to insert, remove and clean; fewer features may fit in the shell.

EARBUD — Hearing-enabled earbuds

Some FDA-authorized products can provide OTC hearing-aid functionality. Confirm compatibility, intended use, battery life and phone requirements.

General style descriptions based on NIDCD and FDA consumer guidance [1][5]. Product capabilities vary.

OTC or prescription?

Both are regulated medical devices, but they serve different users and offer different support models.

	OTC hearing aids	Prescription hearing aids
Intended user	Adults 18+ with perceived mild-to-moderate loss	People of any age and any degree of loss, as indicated
How purchased	In stores or online without a prescription	Through a licensed hearing professional
Fitting	User-controlled; some are self-fitting	Professionally selected, programmed and adjusted
Best fit	Comfortable with technology and self-directed care	Complex needs, severe loss or preference for hands-on care
Support	Varies: app, phone, remote or retail support	Usually includes in-person evaluation and follow-up
Medical exam	Not required for eligible adults, but warning signs still matter	Requirements may vary; medical referral when indicated

There is no universally 'best' route

The right choice depends on symptoms, hearing level, dexterity, comfort with apps, access to care, budget and desired support. Seek professional help whenever the problem is unclear or the first device does not help enough.

Source: FDA OTC and prescription hearing-aid guidance [1][2][3].

WORKSHEET 1

Worksheet: choose your care route

Check the statements that sound like you. Bring this page to a professional if you are unsure.

OTC may be worth exploring

- [ ] I am 18 or older.
- [ ] My hearing difficulty seems mild to moderate.
- [ ] I have none of the warning signs on page 7.
- [ ] I am comfortable following app or device instructions.
- [ ] I am willing to adjust settings and contact support.
- [ ] A lower upfront cost is an important priority.

Professional care may fit better

- [ ] My hearing difficulty is severe, sudden, uneven or unclear.
- [ ] I want a complete hearing evaluation first.
- [ ] I prefer professional fitting and verification.
- [ ] I have dexterity, vision or technology concerns.
- [ ] I need help with tinnitus, dizziness or complex symptoms.
- [ ] Ongoing local support matters more than self-service.

My preferred next step

- [ ] Schedule a hearing evaluation
- [ ] Research OTC options
- [ ] Talk with an ENT
- [ ] Ask my insurer about benefits

Questions I still have

This worksheet supports discussion; it is not a screening or diagnostic tool.

Choosing a hearing aid

Choose features because they solve a real listening problem - not because they sound impressive.

- **Directional microphones**
  Can prioritize sound from the front and may help in some noisy situations.
- **Rechargeability**
  Avoids frequent battery changes, but requires a reliable charging routine and access to power.
- **Telecoil**
  Connects with hearing loops in some theaters, houses of worship, counters and public spaces.
- **App controls**
  May adjust volume, programs, direction and support. Consider privacy, updates and ease of use.
- **Feedback control**
  Helps reduce whistling. Persistent feedback can also signal fit, wax or maintenance problems.
- **Noise reduction**
  May make background sound more comfortable; it does not erase competing voices.
- **Bluetooth streaming**
  Sends calls or media to the aids. Check phone, operating system and app compatibility.
- **Remote microphones**
  Places a microphone near a speaker and can be especially useful across distance or in noise.
- **Rechargeable case**
  A portable charging case can help travelers, but adds size and replacement cost.
- **Water resistance**
  Ratings differ. Water-resistant does not automatically mean safe for showering or swimming.

Shortlist checklist

- [ ] The device fits my top three listening situations.
- [ ] I can insert, remove, charge and clean it comfortably.
- [ ] My phone and required app are compatible.
- [ ] The support, return and warranty terms are clear.

Feature performance and compatibility vary by model. Verify specifications before purchase.

WORKSHEET 2

Worksheet: your listening priorities

Rank the situations that matter most. Use 1 for highest priority, then add notes.

RANK	LISTENING SITUATION	WHAT IS DIFFICULT NOW?
	One-to-one conversation at home	
	Family meals and small groups	
	Restaurants or crowded rooms	
	Telephone or video calls	
	Television and streaming	
	Work meetings or classrooms	
	Music, theater or worship	
	Outdoor activity and exercise	
	Driving or riding in a car	
	Travel, airports and public spaces	

Tip: the best comparison is the one built around your actual day, not a generic feature count.

A smarter buying roadmap

Slow the process down enough to compare clearly, while using trial time efficiently.

BEFORE

Define three listening goals

Use page 14. Be specific: 'follow dinner conversation without asking for repeats' is more useful than 'hear better.'

RESEARCH

Choose a route and shortlist

Compare OTC and prescription support, then shortlist two or three options that fit your needs and handling ability.

VERIFY

Read the package and policies

Confirm intended user, phone needs, batteries, return window, nonrefundable fees, warranty, support and repair process.

BUY

Keep proof and documents

Save the receipt, serial numbers, packaging, user guide, support contacts and the exact date the return period ends.

TEST

Use a structured trial

Try the same real-world situations with each option. Record comfort, clarity, effort, reliability and support quality.

DECIDE

Keep, adjust, exchange or return

Do not wait until the final day. Ask for help early and preserve enough time to complete any return.

Based on FDA purchasing guidance [2][3].

Understand the total cost

The sticker price is only one part of value. Compare what happens after the sale.

DEVICE	Price for one aid or a pair; charging case; required accessories
CARE	Hearing evaluation; fitting; verification; follow-up visits; remote support
UPKEEP	Domes, wax guards, filters, batteries, receivers, cleaning tools
RISK	Return shipping; restocking or professional fees; loss-and-damage coverage
LONG TERM	Repairs after warranty; app/phone compatibility; replacement cycle

Coverage note for U.S. readers

Original Medicare does not cover hearing aids or exams performed solely to fit hearing aids. Medicare Part B may cover diagnostic hearing and balance exams under specified conditions. Medicare Advantage, Medicaid, employer plans, veterans benefits and discount programs vary - confirm details directly with your plan or program before purchase.

Source: Medicare.gov hearing and balance exam coverage [7]. Benefits and rules can change.

My maximum all-in budget: $_____

Read the terms before you pay

A strong return policy and responsive support can be more valuable than an extra feature.

Return or trial period

How many calendar days? When does the clock start? Is advance authorization required? Who pays shipping? Are any fees nonrefundable?

Warranty

What is covered and for how long? Are charger, receivers and accessories covered? Is loss or damage included? Who pays shipping?

Support

Phone, chat, video or in-person? Which days and hours? Is support provided by a licensed professional, technician or customer-service team?

Repairs

Where is the device serviced? Is a loaner available? What are expected turnaround time and out-of-warranty costs?

Compatibility

Which phones and operating systems work today? Is a smartphone required for setup or daily control? What happens if software support ends?

Seller credibility

Is contact information clear? Are claims realistic? Can you obtain full labeling before purchase? Is the return address and process specific?

FDA guidance emphasizes total cost, trial terms, nonrefundable fees, warranty, care and support [3].

Worksheet: compare two options

Fill this out before the trial deadline. Use evidence from your actual listening priorities.

	OPTION A	OPTION B
Brand / model		
OTC or prescription		
Price for pair	$	$
Return deadline		
Nonrefundable fees		
Warranty / repairs		
Support included		
Phone / app fit		
Handling / comfort	1 2 3 4 5	1 2 3 4 5
Speech in quiet	1 2 3 4 5	1 2 3 4 5
Speech in noise	1 2 3 4 5	1 2 3 4 5
Sound quality	1 2 3 4 5	1 2 3 4 5
Support quality	1 2 3 4 5	1 2 3 4 5

Scale: 1 = poor, 3 = acceptable, 5 = excellent. Do not average away a safety or comfort problem.

Your first 30 days

Adjustment is a process. Increase challenge gradually and ask for help early.

DAYS 1-3

Build the routine

- Practice insertion and removal in good light.
- Start in calm, familiar environments.
- Learn charging or battery steps.
- Write down discomfort or confusing controls.

DAYS 4-7

Add everyday sound

- Wear the devices for longer periods as comfortable.
- Notice clarity and effort, not only volume.
- Try one-to-one conversation and television.
- Contact support for pain, poor fit or harsh sound.

WEEK 2

Test your priorities

- Try your top three situations from page 14.
- Record where speech remains difficult.
- Use any programs or directional controls.
- Request adjustments before the return deadline nears.

WEEKS 3-4

Make the decision

- Repeat the same situations after adjustments.
- Review return, warranty and repair terms.
- Confirm reliability, comfort and support quality.
- Keep, exchange or return with time to spare.

FDA advises regular use, practice and testing in varied environments; adaptation takes time [3][4].

Make conversations easier

Hearing aids work best when the environment and communication habits also support you.



CHANGE THE ENVIRONMENT

- Face the person and keep their face well lit.
- Reduce competing television or music when possible.
- Move closer rather than only increasing volume.
- Choose a seat with your back toward a wall in restaurants.

CHANGE THE CONVERSATION

- Ask for a rephrase, not just a louder repeat.
- Tell people which side or distance works best.
- Use captions for television and video calls.
- Consider a remote microphone for distance or group situations.

A useful phrase

'I want to hear you. Please face me and say that another way.'

Image is illustrative. Communication strategies can complement, but do not replace, appropriate hearing care.

Daily care protects performance

Follow your model's instructions. Small maintenance habits prevent many common problems.



EVERY DAY

Wipe the device with a clean, dry cloth. Check openings for visible debris. Charge or replace batteries as instructed.

REGULARLY

Replace domes, filters or wax guards on the manufacturer schedule or when performance changes.

KEEP DRY

Protect from showers, swimming, high humidity and extreme heat unless the manufacturer explicitly allows exposure.

STORE SAFELY

Keep aids, batteries and small parts away from children and pets. Use the case when not wearing them.

AVOID

Do not use water, solvents or alcohol on the electronics unless the instructions specifically permit it.

INSPECT

Look for cracks, blocked tubing, damaged wires, poor charger contact or loose domes. Contact support when unsure.

Source: FDA care guidance [3]. Always follow the instructions for your specific device. Image is illustrative.

Troubleshooting: start simple

If a problem is painful, sudden or persistent, stop and contact qualified support.

PROBLEM	CHECK FIRST	NEXT STEP
No sound	Power, charge/battery, correct left/right aid, blocked opening	Restart; clean only as directed; contact support
Sound is weak	Wax guard/filter, dome, battery/charge, program and volume	Replace approved parts; request adjustment
Whistling	Insertion, dome/earmold fit, wax, nearby hat or hand	Reinsert; contact a professional if persistent
Harsh or too loud	Program, volume, fit and recent settings changes	Reduce as instructed; request prompt fine-tuning
Bluetooth drops	Phone compatibility, Bluetooth state, app permissions, distance	Restart and re-pair using official instructions
Ear is sore	Placement, dome size, skin irritation or trapped moisture	Stop use if painful; contact professional or manufacturer
Still cannot follow speech	Real-world fit, device power, settings, degree/type of loss	Seek professional evaluation; do not only raise volume

Report injuries, overheating, lodged pieces or sudden worsening associated with a device to a clinician and the FDA MedWatch program [3].

For family and communication partners

Good communication is shared work. Support should protect dignity and independence.

Helpful habits

- Get the person's attention before speaking.
- Face them and keep your mouth visible.
- Speak clearly at a natural pace.
- Rephrase if repetition does not help.
- Reduce distance and background noise.
- Include the person; do not answer for them.
- Confirm important details in writing.
- Ask what works instead of assuming.

Avoid these common mistakes

- Shouting from another room.
- Speaking while turning away.
- Saying 'never mind' after a missed phrase.
- Testing the person unexpectedly.
- Assuming the device should create perfect hearing.
- Treating support needs as a lack of intelligence.

Plan together

Choose one shared goal for the next two weeks - for example, a more comfortable family dinner - then change one factor at a time and notice what helps.

WORKSHEET 4

Worksheet: your hearing action plan

Turn what you learned into a short, specific next step.

My three most important listening goals

Warning signs or medical questions to discuss

My preferred care route and why

My maximum all-in budget

Two options I want to compare

My next action, owner and date

Keep this page with your receipt, test results and trial deadline.

Plain-language glossary

Useful terms you may see in product pages, packages or appointments.

Audiogram

A graph of the quietest tones heard at different pitches.

Audiologist

A hearing professional trained to evaluate hearing and provide non-medical treatment, rehabilitation and hearing technology services.

ENT

An ear-nose-throat physician, also called an otolaryngologist, who diagnoses and treats medical or surgical ear conditions.

Decibel (dB)

A unit used to describe sound level. Audiograms often use dB Hearing Level (dB HL).

Frequency (Hz)

Pitch. Lower numbers are lower sounds; higher numbers are higher sounds.

Feedback

Whistling that can occur when amplified sound re-enters the microphone, sometimes because of fit or blockage.

OTC hearing aid

An FDA-regulated hearing aid for adults 18+ with perceived mild-to-moderate hearing loss, available without a prescription.

Prescription hearing aid

A hearing aid obtained through a licensed professional and appropriate for broader ages and degrees of loss, as indicated.

PSAP

Personal sound amplification product for people with normal hearing in specific situations; not a hearing aid.

Real-ear measurement

A professional verification method that measures sound near the eardrum while the hearing aid is worn.

Receiver

The miniature speaker that delivers amplified sound into the ear.

Telecoil

A small coil that receives audio from compatible phones or hearing-loop systems.

Tinnitus

Perception of ringing, buzzing or another sound without an external source.

Wax guard

A replaceable filter that helps keep earwax out of a receiver or sound outlet.

Definitions are simplified for consumer education. Ask your professional how terms apply to your situation.

Sources and further reading

Primary consumer sources accessed August 25, 2026.

1 U.S. Food and Drug Administration
Hearing Aids
https://www.fda.gov/medical-devices/consumer-products/hearing-aids

2 U.S. Food and Drug Administration
OTC Hearing Aids: What You Should Know
https://www.fda.gov/medical-devices/hearing-aids/otc-hearing-aids-what-you-should-know

3 U.S. Food and Drug Administration
How to Get Hearing Aids
https://www.fda.gov/medical-devices/hearing-aids/how-get-hearing-aids

4 U.S. Food and Drug Administration
Hearing Aid Benefits and Limitations
https://www.fda.gov/medical-devices/hearing-aids/hearing-aid-benefits-and-limitations

5 National Institute on Deafness and Other Communication Disorders
Hearing Aids
https://www.nidcd.nih.gov/health/hearing-aids

6 National Institute on Deafness and Other Communication Disorders
Age-Related Hearing Loss
https://www.nidcd.nih.gov/health/age-related-hearing-loss

7 Medicare.gov
Hearing and Balance Exams
https://www.medicare.gov/coverage/hearing-balance-exams

8 HearingAids.md
Compare Hearing Aids
https://hearingaids.md/compare

9 HearingAids.md
Hearing-Aid Brand Directory
https://hearingaids.md/brands

10 HearingAids.md
How to Read Your Audiogram
https://hearingaids.md/blog/understanding-audiogram-results

Editorial policy

Health and regulatory guidance can change. Before each new edition, verify every time-sensitive statement, link, product description and coverage note against current primary sources. Do not display a medical-review badge unless a qualified reviewer has completed and documented that review.

Notes

Use this space for appointment questions, support contacts, serial numbers and trial dates.

PROFESSIONAL / CLINIC

MANUFACTURER SUPPORT

PURCHASE DATE

RETURN DEADLINE

WARRANTY END DATE

SERIAL NUMBERS

Questions and observations

HearingAids.md



Learn. Compare. Understand. Choose.

Independent hearing-health education, practical comparisons and clear consumer guidance.

VISIT HEARINGAIDS.MD

Version 1.0 | August 2026 | Copyright 2026 Sonivue Lab

For educational purposes only. Not medical advice.